

LAB-K

LABORATORY CERAMIC WORKTOP



Corrosion
Resistant



Scratch
Resistant



Fireproof
up to 1200 °C



Antibacterial




Non-deformation



Easy to clean



Non-inflammable
Class 1A



Waterproof



White



Grey



Black

Available in 2600 x 800 x 15mm

Lab-K laboratory 15mm ceramic worktop are sintered at a high temperature of 1200 °C can truly achieve the same color of the body and glaze and avoid problems such as delamination and shedding. Ceramic worktop is corrosion-resistant, fireproof, waterproof, scratch-resistant, abrasion-resistant, non-deformation.

Ceramic worktops can resist strong acid, alkaline and organic solvent corrosion (except for hydrofluoric acid equivalent chemical reagents).



Laboratory work bench with Black and Grey colour ceramic worktop



KL LAB SUPPLIES SDN. BHD. (718086-W)

No. 5, Jalan 36/41, Off Jalan Dato M. Sellathevan, Batu 4, Seksyen 36, 40470 Shah Alam, Selangor, MALAYSIA.

Tel : +603 5167 3601 +6011-3585 1691 /2 /3 Fax : +603-5167 3602 H/P : +6012-202 7055

Website : www.kllab.com Email : info@kllab.com