

Lab-K

Polypropylene Acid Safety Storage Cabinet



- ❖ **Lab-K** Acid Safety Storage Cabinet is all 10mm thick polypropylene constructed.
- ❖ Polypropylene is the material of choice for acid / corrosive chemical storage for laboratory environments and particle free cleanroom storage.
- ❖ You will enjoy years of clean, safe storage of your hazardous acid / corrosive chemicals and materials without the rusting due to metal construction or the delaminating of wood core products.

Model : CA212-451

Size : 900mm (W) x 450mm (D) x 1800mm (H)



Model : CA213-451 Size : 900mm (W) x 450mm (D) x 1800mm (H)

Code No.	Door Type	Capacity	Dimension	Shelves
CA212-451	2 Door Manual Self-Close	210 Liters / 45 Gallons	W = 900 mm D = 450 mm H = 1800 mm	4
CA213-451	4 Door Manual Self-Close	210 Liters / 45 Gallons	W = 900 mm D = 450 mm H = 1800 mm	3


Product Features



- ❖ Door with glass panel for easy selection. This will reduce the acid fume escape due to door opening.
- ❖ Cabinets are lockable to prevent unauthorized access.



- ❖ Adjustable shelves with vertical lip for safety.
- ❖ Spillage tray to base area to capture incidental leaks. These trays can be easily removed and cleaned as required.
- ❖ Designs with Exhaust venting on top of the cabinet.



KL Lab Supplies Sdn. Bhd. (718086-W)

No.5, Jalan 36/41, Off Jalan Dato M. Sellathevan, Batu 4, Seksyen 36, 40470 Shah Alam, Selangor, MALAYSIA

Tel : +603-5167 3601 / +6011-3585 1691/ 2/ 3 H/P : +6012-202 7055

Website : www.kllab.com Email : info@kllab.com